



A TOOL FOR NANO-MILLING POORLY WATER-SOLUBLE API USING THE PRECELLYS® + CRYOLYS

Bertin Instruments, Montigny-le-Bretonneux

CONTEXT

The majority of approved drugs have poor solubility and permeability. Reducing the particle size and surface area of the API (Active Pharmaceutical Ingredient) **increases drug performance and bioavailability.**

The Precellys® 24* with Cryolys was validated for drug nano-sizing using different milling protocols.

MATERIALS & PROTOCOLS

Materials

- Precellys® 24* Homogenizer + Cryolys
- Precellys lysing kits: different bead sizes, material, and concentration were tested in 2mL reinforced tubes
- Samples: 5 mg/mL of Indomethacin, Ibuprofen, and Itraconazole were used as models of poorly watersoluble API
- Aqueous vehicles: 0.5% Pluronic® F68 and 0.5% HPC (low viscosity) dissolved in 1mL of ultrapure water
- Oily-based vehicles: 1mL each of soybean and peanut oils (Olvea), glycerol monocaprylocaprate (Capmul® MCM) and triglycerides (Miglyol® 812N and 840)



Protocols

- Precellys® 24*: 5500 rpm, 24 cycles x 30s, 5s pause
- Cryolys cooling unit: Maintained a sample chamber temperature between 2-4°C
- Nano-milling efficiency was studied using different parameters: bead diameter (0.5-2.8 mm), bead concentration (1-4 g/ml) and material (yttrium stabilized zirconia, zirconium oxide, VHD zirconium oxide and stainless steel). Particle size after nano-milling was measured by a Malvern Zetasizer (dynamic light scattering).

* technically equivalent to Precellys 24 Touch



A TOOL FOR NANO-MILLING POORLY WATER-SOLUBLE API USING THE PRECELLYS® + CRYOLYS

RESULTS

The particle size distribution of indomethacin (IN), ibuprofen (IB), and itraconazole (IT) was determined using **microscopy**. The average particle diameter (D50) was calculated as 27 µm (IN), 144 µm (IB) and 14 µm (IT), respectively. When using aqueous vehicles, the smallest particle sizes were obtained using 0.6- 0.8mm yttrium stabilized zirconia at a bead concentration of 2 grams/mL (Figure 1A, highlighted in yellow). The use of Capmul MCM® yielded the smallest particle sizes for all three drug compounds tested (Figure 1B, highlighted in green).

API	IN	IB	IT	IN	IB	IT	IN	IB	IT
Aqueous vehicle (1mL)	0.5% Pluronic F88 and 0.5% HPC								
Bead Concentration and Type	2 g/mL of 0.5mm VHD ZrO			2 g/mL of 0.6-0.8mm yttrium stabilized zirconia			2 g/mL of 1-1.2mm ZrO2		
D ₅₀ (nm)	380.6	677.1	1426	326.8	464.3	719.2	348	601	700

Oily vehicle (1mL)	Soybean oil			Peanut oil			Capmul® MCM			Mglyo® 812N			Mglyo® 840		
API	IN	IB	IT	IN	IB	IT	IN	IB	IT	IN	IB	IT	IN	IB	IT
D ₅₀ (nm)	95 + (2967 - 38%)*	456	497	500	433	815 + (5172 - 1%*)	150 + (4737 - 11%*)	275	373	489	409	1024	633	534	1767 + (5412 - 2%*)


Figure 1A. Poorly water soluble API were reduced to nano-size using small zirconia beads in aqueous stabilizers. The best results were obtained with bead sizes ranging between 0.5 – 1.2mm at a bead concentration of 2 grams/mL. IN=Indomethacin; IB=Ibuprofen; IT=Itraconazole. Particle size in red indicates a non-homogeneous mixture.

Figure 1B. Poorly water soluble API were reduced to nano-size using 0.6-0.8mm yttrium stabilized zirconia at a bead concentration of 2 grams/mL in 5 different oily-based vehicles. *represents the % intensity of larger particle sizes in the sample when 2 peaks are observed.

CONCLUSION

Nano-milling of poorly water-soluble drugs can be easily performed using the Precellys® 24* + Cryolys in both aqueous and oily based vehicles, while limiting API consumption. Amongst the different dispersing media tested, Capmul® MCM yielded the smallest particle size, and reduced indomethacin to 150 nm. The Precellys® + Cryolys represents a new technology for drug nano-milling that can be utilized in early discovery and preclinical research.

DISCOVER OUR COMPREHENSIVE RANGE OF SOLUTIONS



Join the Bertin Technologies community!

1500+ documents available online on the Application center

www.bertin-technologies.com/application-center/

* technically equivalent to Precellys 24 Touch