

Precellys Connect



Comprehensive mobile application for Precellys Evolution Touch

- Live monitoring of the instrument parameters
- Direct access and download of protocols from Bertin Technologies database
- Protocol creation and transfer directly on the App
- Personal database protocol management on the App
- Available in English & French

PRECELLYS CONNECT

Precellys Connect is a dynamic mobile app designed to remotely monitor and program your Precellys Evolution Touch instrument. Those functionalities are reachable once the bluetooth connection is established between the instrument and your device. The App allows you searching protocols from the extensive Bertin Technologies database and then to transfer them directly or once modified to your Precellys Evolution Touch.

Optimize your workflow by integrating your protocols directly from your smartphone into your instruments and vice versa, improving productivity and efficiency.

Incorporate Precellys Connect into your daily routine and R&D projects to optimize the performance of your Precellys Evolution Touch instrument. Get the full potential of your experiments thanks to the flexibility offered by Precellys Connect.

Quick User Guide



1. Download & Install the application from the website or scan the QR code



2. Connect your device with the desired instrument using the application



3. Set and load your experiment on the instrument



4. Monitor the process remotely

Precellys Connect is available on the Precellys Evolution Touch (contact us if you want to check the compatibility of your instrument).
Precellys Connect is only available on Android devices.

Discover our comprehensive range of solutions



AIR
SAMPLING



SAMPLE
PREPARATION



ASSAY
KITS



TECHNOLOGY
ON-DEMAND



Bertin Technologies SAS
Parc d'activités du Pas du Lac
10 bis, avenue Ampère
78180 Montigny-le-Bretonneux, FRANCE
01 39 30 60 00

WWW.BERTIN-TECHNOLOGIES.COM