


Multi Tissue Dissociation Kit



The reliable kit for obtaining single viable cells

- Bead-shearing dissociation
- High viability rate
- Versatility

GET HIGH RATE OF LIVING SINGLE CELLS FROM YOUR STARTING TISSUE

Our Multi Tissue Dissociation Kit combines enzymatic and mechanical dissociation.

This approach maximizes cell viability and significantly increases the yield of living single cells from various tissue types.

The kit includes a specially formulated **enzyme mix** and **15 mL tubes**, for up to **20 preparations**.

Key Features:

- > 80 % viability rate
- Patented Bead-Shearing Technology
- Versatility



D05717
Multi Tissue Dissociation Kit




D30717
Dissociation Tube

Validation on a Range of Tissues for Reliable Results

The Multi Tissue Dissociation Kit has been validated on tissues from:



Liver
(Mouse)



Spleen
(Mouse)



Lung
(Mouse/Human)



Thymus
(Mouse)



Heart
(Mouse)

This kit is meant to be use with Bertin Technologies Precellys Evolution instrument

VISIT OUR WEBSITE



Bertin Technologies SAS
Parc d'activités du Pas du Lac
10 bis, avenue Ampère

78180 Montigny-le-Bretonneux, FRANCE
01 39 30 60 00

EU webstore: bertin-bioreagent.com | US webstore: bertin-corp.com

WWW.BERTIN-TECHNOLOGIES.COM